


## STEP 1 - Turn ActiPatch ON

### To turn on ActiPatch.

Grasp the patch firmly as shown and pull out the white key protruding from the patch.




The red indicator light confirms that the device is operating properly. To avoid sleep disruption, the brightness of the red indicator light is minimized and is only visible in a dark location.

This device is designed to deliver 720 hours of effective therapy.

## STEP 2 - Place ActiPatch:

Place ActiPatch directly on the injured tissue, or on top of the bandage or wound dressing covering the injury.

- Place the Actipatch so that the area for treatment is located directly below and inside of the antenna perimeter, as indicated by blue dots in the illustration.



- The patch should be placed so that the flat side is toward the skin.

## STEP 3 - Affix ActiPatch:

### Using Medical Grade Tapes

- Place ActiPatch directly over the injured tissue.
- Tape ActiPatch in place to hold the antenna directly over the injured tissue. You may use cloth, paper, transparent or wound closure tape.
- Remove tape gently to avoid tearing the foam covering.

### Using Elastic Bandages or Cohesive Wraps

- Place ActiPatch over injured tissue.
- Wrap ActiPatch in place to hold in the proper position.

Caution: DO NOT OVER-TIGHTEN THE ELASTIC BANDAGE OR COHESIVE WRAP AS THIS MAY RESTRICT CIRCULATION.

The effectiveness of ActiPatch is affected by the elapsed time between the occurrence of trauma and first use of the unit and the initial duration and consistency of follow on application. If possible apply 24 hours per day or minimally 6 to 8 hours daily.